

SPIRARE®



Emergency

- Ships
- Offshore
- Platforms
- Ambulances

Module Emergency ECG

The Spirare Emergency ECG module is an optional module for Spirare. This module is developed for emergency ECG recordings on patients on ships, offshore platforms or in mobile units such as ambulances.

Using an ECG belt with chest electrodes enables correct connection of patient within 20 seconds. Standard ECG electrodes and cables are also available.

The Spirare software is straight forward to use and requires minimal training to operate.

ECG recordings are sent by email (as attachments) to a cooperating hospital or medical unit for evaluation. Patient information as well as images of the patient can easily be included in the email.

The recipients in a hospital or another medical unit, using Spirare, may store the transferred ECG data in a database.

System requirements: Microsoft Windows XP/2000/2003/Vista and a working email program. The email addresses of two recipients may be preregistered in Spirare.

DIAGNOSTICA AS
Sondreveien 1
N-0378 Oslo
Norway

Tel (+47) 22 92 40 00
Fax (+47) 22 92 40 09

diag@diag.no

the standard in network ECG