



Portable Unit

- Ambulatory use
- Database replication

Portable Unit Module

The Portable Unit module makes Spirare the complete ambulatory system for spirometry and ECG.

Perfect for

- Occupational health screening
- Bedside testing in hospitals
- Ambulatory use by general practitioners or nurses
- Rehabilitation home care units

The Portable Unit module allows transfer of patient records from the main office database to a portable PC. The PC can subsequently be used separately (offline) to perform tests that later may be replicated back to the main office database, making the data immediately accessible for all Spirare users.

The module is scalable and permits several portable units to be used with one database.

Portable unit options

- Synchronise data on the portable unit with the main office database
- The portable unit database is easily removed.

Both manual and scheduled database synchronisation are possible, to accommodate a wide range of network installations.

– Spirare necesse est!

the standard in network spirometry and ECG

DIAGNOSTICA
Sondreveien 1
N-0378 Oslo
Norway

Tel +47 2292 4000
Fax +47 2292 4009

diag@diag.no